


### Center-open Vise Grip / GTW-50

A sample, regardless of the size, is always gripped in the center of the vise as the grips open and shut at the center.




Lever type

[Features]

- The grips always hold a test sample at the center of the vise.
- The jaw faces have a vertical V-groove that can firmly hold a round shaft.
- Easy-to-operate dial knob type and strong-grip lever type are available.

[Dimensions]



Dial knob type

[Specifications]

Type	Capacity	Jaw face width	Weight
GTW-50R	500N	50mm	3kg
GTW-50L	1000N	50mm	3kg

[Note]

- Some samples may not correctly be gripped depending on material and shape.
- Information within this document is subject to change without notice.